

# Burnaid®

## BURNAID PRODUCTS

### What does Burnaid

- Relieves pain.
- Cools by pulling heat out of the burn.
- Sterile barrier against contamination.
- Non-adhesive.
- Antimicrobial - helps prevent infection.
- Cooling relief up to 2 hours.
- Ideal for transport or where larger burns may occur.
- Compresses to cover all parts of the body - burns, hands, limbs, trunk, face.

*NOTE Burns greater than a two pence coin or deep should always exhibit for doctors.*

### Burnaid - Gel

First aid for burns, sunburn, scalds, insect bites and minor injuries. In all places where a minor burn may occur.

#### Burnaid spray 50ml

- Spray liberally over the injury.
- Not to be rubbed!
- Allow the gel to be on for 30 - 50 min.

#### Burnaid gel 3,5g

- Tear off and squeeze out the gel on the burn.
- Not to be rubbed!
- Allow the gel to be on for 30 - 50 mins.
- If damage to the fingertips, insert your finger into the bag with gel..



#### Dressing sizes:

10x10, 20x20, 30x40 (face mask), 55x40, 120x200 (body) cm

### Burnaid – Dressings

First aid for burns and scalds  
At all locations where major burns can occur.

#### Burnaid dressing

- Tear off and remove the dressing and place on the burn.
- Place the empty bag over the dressing and gently wrap a bandage over to keep the dressing in place.
- Allow the gel to remain on for up to 2 hours.