




PRODUKCATALOGUE - 2023



AlphraMedic AB - 19 year
2004-2023



AlphraMedic AB has extensive experience from sales and manufacturing of products to different users. AlphraMedic manufactures and supplies products within First Aid to retailers in Scandinavia and the rest of Europe.

AlphraMedic AB manufactures i.e., different types of holders, dispensers and kits for First Aid products. The production consists of partly in-house developed products and partly products manufactured according to customer wishes.

AlphraMedic AB's business concept is to satisfy our customers with a wide range of First Aid products such as:

Dispensers and Holders, Plasters on a roll or Plaster strips, Gel dressings, Burn products, Blood stop spray, Wound plaster spray, First aid bags (soft/hard), Refill range, etc.

Contact us if you have products that must be customized to best suit the market and your customers!

Contact us if you would like more information about products and retail prices!

With Best Regards
AlphraMedic AB



Addr: Kämpevägen 17, S-553 02 Jönköping,
Sweden

Phone: + 46 36 340485

Mob: +46 70 273 77 40

Mob: +46 70 993 99 10

E-mail: info@alphramedic.se

Web: www.alphramedic.se



PLASTERDISPENSER

to 6x500 cm roll


Suites all types of 5m plaster roll

The D-56 is a popular Dispenser for easily cutting the size of patch you need, wide or narrow, you choose!

A practical, simple, economical and flexible system! Mounted on the wall!

Made of clear acrylic plastic, width 6cm.

Ref. No. 3000-1 D-56 dispenser



Suitable for all kinds of plaster rolls, e.g. Textile, Non Woven, Polyurethane foam, etc

PLASTER ON ROLL – 5 METER


Wide: 6cm and 8cm



- **Textile patch - 6cm**
- Use on clean and dry skin
- Elastic and comfortable to use
- Painless removal of the patch
- Fits D-56 plaster dispenser
- Cut to desired size

Ref. No. 3361

- **Blue Detectable Textile with metal foil - 6cm**
- Use on clean and dry skin
- The blue color provides visual identification
- The wound pad with metal foil can be identified with a metal detector
- Elastic and comfortable to use
- Painless removal of the patch
- Fits D-56 plaster dispenser
- Cut to desired size



Ref. No. 3360



- **A strong textile patch - 8cm**
- Use on clean and dry skin
- Elastic and comfortable to use
- Painless removal of the patch
- Cut to desired size (cut off)

Ref. No. 3361-8

- **Non-Woven patch - 6cm**
- Hypoallergenic
- Suitable for sensitive skin
- Use on clean and dry skin
- Elastic and comfortable to use
- Painless removal of the patch
- Fits D-56 plaster dispenser
- Cut to desired size

Ref. No. 3362



SUPERSOFT


Wide: 6cm

Length: 450cm

SuperSoft - Blue


- Self-vulcanizing foam plastic bandage
- The blue color provides visual identification
- Elastic and comfortable to use
- Does not stick to wounds or hair
- The bandage breathes, which provides good healing
- Fits D-56 plaster dispenser
- Cut to desired size
- For fixation of burn compresses

Ref. No. 3385-B (blue)




Ref. No. 3000-1

1



2



3




NEW

EASYON DISPENSER

Plaster dispenser with patch refill



Ref. No. 2505 - complete



Ref. No. 2510- empty


New EasyOn patches are easy and convenient to use!!

EasyOn dispenser is made of transparent and durable plastic. Used with EasyOn patch refill. See below. With a simple hand grip, you change the refill. Mounted on the wall.

The EasyOn refill fits the EasyOn dispenser and you can change the refill with a simple hand grip. Gently pull down on the bottom edge of the patch and the upper part of the paper on the back of the patch will come off - ready to mount!

EasyOn refill is available in two different materials and three different sizes

Blue Detectable Textile and Beige Textile



25x72mm,

56x70mm,

Fingertip (65x40mm)

Ref. No.	Product	Size	Patch/refill
3360-25-1R	Blue Detec Textile	25x72mm	36
3360-52-1R	" "	56x70mm	18
3360-FT-1R	" "	Fingertip	18
3361-25-1R	BeigeTextile	25x72mm	36
3361-52-1R	" "	56x70mm	18
3361-FT-1R	" "	Fingertip	18



25x72mm

25x72mm

65x40mm Fingertip

*Delivered in a box with 5 or 6 refills (alt. 10 refills) Not bulk!
The possibility of own patches is available for larger quantities!*



FIRST AID CASE – M **NEW**

First aid bag with handle - red



Size: 21,5 x 7,5 x 17,5 cm

A stylish first aid bag with a handle and an inner pocket that you can fold out. Place for an ex. DIN kit, or if you want to add Burncare products. Or other products such as Blood Stop - o Wound plaster spray, Pharmagel, Suture tape, Surgical tape, Waterproof plaster, Plaster refill, Long plaster, Vendapress - wrap, Wound tissue, etc.

Ref. No. 5558



FIRST AID CASE – XL **NEW**

First aid bag with handle and shoulder strap - green



Size: 27,5 x 25,5 x 12,5 cm

A large first aid bag with handle and shoulder strap. Place for an ex. DIN kit in one half and in the other half e.g. Burncare products and other products such as Blood Stop - o Wound Plaster spray, Eyewash 250ml alt. 50ml, Pharmagel, Suture tape, Surgical tape, Waterproof plaster, Plaster refill, Long plaster, Winding press - wrap, Wound wipe, etc.

Ref. No. 5733



POCKET BREATHING MASK with valve

Practical pocket packaging with cabin hook for flexible attachment

Protection against blood contamination and infection during the mouth-to-mouth method

Instructions for use

1. Clear any obstructions in the patient's airways
2. Tilt the patient's head back
3. Unfold the breathing mask and place it with the green oval valve in the patient's mouth
4. Make it as tight as possible between the mask and the patient's face
5. Start the word of mouth method

Contains: 1 breathing mask with valve
2 protective gloves

Ref. No. 6053



BREATHING MASK with valve

Protection against contamination during the mouth-to-mouth method

Instructions for use

1. Clear any obstructions in the patient's airways
2. Tilt the patient's head back
3. Unfold the breathing mask and place it with the green oval valve in the patient's mouth
4. Make it as tight as possible between the mask and the patient's face
5. Start the word of mouth method



Ref. No. 6010



WOUND PLASTER SPRAY

Wound dressing spray is an aerosol with quick-drying wound dressing that provides quick and safe protection. It is made to protect minor wounds and injuries such as lacerations, abrasions and cuts at workplaces, playing sports and at home.

Wound plastic spray reduces the risk of infection, protects against water and dirt, which provides safe showering and staying in wet environments or dusty and dirty workplaces, simple, convenient and safe to use, good in areas where it is difficult to use ordinary plasters and bandages.

Ref. No. 3140 60ml



BLOOD STOP SPRAY alginat

BLOOD STOP-alginate spray is an aerosol with blood-stopping powder. The active substance - calcium alginate - forms a powder layer over the injury. The alginate layer absorbs blood from the wound and exudate forms, forming a cooling, gel-like mass, which aids in the healing process. Chamomile and calendula extract Seaberry oil, Aloe Vera and Jojoba wax provide a natural moisturizing of the damaged area and help in the healing process.

It is made to stop minor bleeding from superficial injuries such as lacerations, abrasions and cuts at workplaces, playing sports and at home.

Press the spray button for a few seconds and cover the damaged area. Spray from a distance of approx. 15-20 cm. For repeated treatment, shake the bottle thoroughly each time.

Ref. No. 3130 60ml



Burn plaster

Pharmagel hydrogel dressing - moist wound healing

Treatment of burns and after-treatment of burn injury*

- Pharmagel is a transparent polyurethane gel compress - 60% water
- Pharmagel Burn plasters cool, relieve and reduce pain
- **ATTENTION!** Wash the wound clean before application. It must be Dry and Clean!
- The healing process can be followed without removing the dressing
- The compress has a self-adhesive polyurethane film
- Apply to the burn, covering the entire wound surface
- The compress can be applied up to 1 week
- Moist wound healing accelerates the formation of new tissue
- Minimizes scarring and scabbing
- Does not stick to wounds and has a good dampening effect
- Individually wrapped and sterilized
- Allows free passage of gases and water vapor
- Protects the skin from dehydration
- Pharmagel is waterproof, enabling safe showering



* A **burn** is primarily treated with **Burncare** gel and gel compresses for up to **two** hours. (alt. cool water)

With **Burncare** you achieve three advantages:

- 1) Cooling to minimize tissue damage,
- 2) Pain relief for increased comfort
- 3) Reduces the risk of additional offensive infection.

After this, **Burncare** has "done its thing". For continued wound healing, **Pharmagel** compresses are used

Smaller burns can of course be used directly with Pharmagel for quick treatment and relief



Size 45x65mm; 75x100mm
Ref. No. 4565-C; 4575-C

STERILE EO

CE
0459

AlphaMedic

Eyewash spray

EyeWash - trippel

Buffered sodium chloride solution 9% **NEW**



EyeWash spray eye-shower for the treatment of eye damage, burns and wound washing.

The spray contains 250ml of buffered sodium chloride solution 9mg/litre.

Flushing time - 15 minutes / spray



Ref. No. 8700

Eye damage:

In case of corrosive eye injuries, it is very important to act quickly treatment, preferably within a few seconds.

Flushing time for acid corrosion damage is normally about 5 minutes. Flushing time for alkaline corrosion damage is at least 20-30 min.

In the case of serious corrosion damage, it is important to quickly seek medical attention and then it is important to rinse all the way during transport to the eye-Doctor. It is therefore necessary to include sufficient flushing time.

AlphaMedic's **EyeWash triple** is a bag that contains 3 eye shower sprays that provide 45 minutes flushing time. It is easy to quickly remove the bag from the wall bracket and take it with you during transport.

At the workplace, the case must be within "**reach**" where there is a risk of that an eye injury may occur. It must be positioned so that you can reach it within a couple of seconds.

ATTENTION! When the spray is opened and the seal is broken, the spray should only be used for washing wounds or cooling burns. As the seal is broken and the spray used, full flushing time cannot be guaranteed.



EyeWash-3 + DIN kit 13164



Ref. No. 8710

Ref. No. 8700 Complete case – incl. wall bracket
Ref. No. 8700-E Without Eyewash spray, " " "

AlphaMedic

STERILE R

CE 0413

Eyewash spray

EyeWash Spray 250ml

Buffered sodium chloride solution 9%

EyeWash spray eye-shower for the treatment of eye damage, burns and wound washing.

The spray contains 250ml of buffered sodium chloride solution 9mg/litre.

Flushing time - 15 minutes / spray

Single holder



Triple holder



Eye injury:

In case of corrosive eye injuries, it is very important with fast treatment, preferably within a few seconds.

Flushing time for acid-corrosive damage is normally about 5 minutes
Flushing time for alkaline burns is at least 20-30 minutes.

In the event of severe burns, it is important to seek medical attention quickly and then it is important to flush all the way during transport to the eye doctor. It is therefore necessary to include enough flushing time.

EyeWash-Spray is suitable for use with wall mounts.

Available as Single Holder or Triple Holder. This means that it is Easy to quickly remove the spray from the holder and take it with you damaged and in transit.

In the workplace, eye wash should be within "reach" where there is a risk of that an eye injury may occur. It shall be positioned in such a way as to can reach it within a couple of seconds.

NB! When the spray is opened and the seal is broken, the spray should only be used for Wound Wash or Cooling of Burn. Since the seal is broken and the spray is used, full flushing time cannot be guaranteed.



Ref. No. 8740 Eyewash spray 250ml – refill

Ref. No. 2050 Single holder 52mm

Ref. No. 2060 Triple holder 52mm



STERILE R

CE
0413



Eyewash spray

EyeWash Spray 50ml

Buffrad Natriumklorid lösning 9%

EyeWash spray eye-shower for the treatment of eye damage, burns and wound washing.

The spray contains 250ml of buffered sodium chloride solution 9mg/litre.

Flushing time - 15 minutes / spray

Single holder
35mm



Eye injury:

In case of corrosive eye injuries, it is very important with fast treatment, preferably within a few seconds.

**Flushing time for acid-corrosive damage is normally about 5 minutes.
Flushing time for alkaline burns is at least 20-30 minutes.**

In the event of severe burns, it is important to seek medical attention quickly. And then it is important to flush all the way during transport to the eye doctor. It is therefore necessary to include enough flushing time.

EyeWash-Spray is suitable for use with wall mounts.

Available as Single Holder. This means that it is easy to quickly remove the spray from the holder and take it with you to the injured person and during transport.

In the workplace, eye wash should be within "reach" where there is a risk of that an eye injury may occur. It shall be positioned in such a way as to can reach it within a couple of seconds.

NB! When the spray is opened and the seal is broken, the spray should only be used for Wound Wash or Cooling of Burn. Since the seal is broken and the spray is used, full flushing time cannot be guaranteed.

Ref. No. 8745 Eyewash spray 50ml – refill

Ref. No. 2040 Single holder 35mm



STERILE R



0413



BAGG FIRST AID KIT



- Original Design / Registered Trademark
- Durable, lightweight 420 denier polyester.
- Weather- and water-repellent material.
- Extra attachment points (Daisy chains).
- Double-sided "Heavy Duty" zipper for easy opening and closing.
- Sturdy plastic handle with extra grip.
- Reinforced belt with attachment point.
- Extra "Molle" attachment point.
- Interior divided into several sections to keep contents in the right place.
- 2 large closable sections with zippers and 2 elastic sections.
- Dimensions: 21 x 15 x 11cm.

BAGG first aid kits are provided filled or empty. The content can be easily adapted to your customers' wishes.

Ref. No. 5510– Orange (empty)
Ref. No. 5515– Orange (complete with DIN
kit 13164)

EC - FIRST AID KITS

Lacerations - Burns - Eye injuries

Stylish and practical quality bags that are delivered complete, custom-made or empty. Perfect at the workplace where you can build a complete functioning system for first aid treatment. The bags can be mounted next to each other or underneath each other.

When mounting on the wall, the entire bag is easily accessible to the injured person. Choose shelves for clarity. Of course, it can also be used without a wall mount when needed.

The bags can also be used for other products such as protective gloves, goggles, respirator mask, hearing protection, etc.

SIZES

Large; Medium; Small

DIMENSIONS:

Large: 30x33x12
 Medium: 30x23x12
 Small: 15x23x8

ACCESSORY:

Bracket L / M, S
 Shelves L, M

COLOR:

Green, Blue, Red (Burn Kit),
 Transparent (Eyewash Kit)

Transparent - Empty

Ref. No. 7230 Medium
 7240 Large




Shelves

Ref. No. 7410 Medium
 7400 Large



Ref. No. - Empty

7200+color Small
 7210+color Medium
 7220+color Large



Bracket

Ref. No. 7300 Small
 7310 Large o. Medium,



EC - FIRST AID KITS

DIN+EASYON

Large

CONTENT:

DIN KIT 13164 + EasyOn dispenser with two plasters refill incl. 2 shelves

Accessories: Wall bracket



DIN KIT

EasyOn dispenser
1 Textil-refill + 1 Blue
detectable Textil refill

Table of contents

DIN 13164

- 1 piece of Tape 2,5cmx5m
- 14 Plaster assortment
- 1 First Aid Bandage 6x8cm
- 2 First Aid Bandages 8x10cm
- 1 First Aid Bandage 10x12cm
- 1 Burn compress 60x80cm
- 6 Compresses 10x10cm
- 1 Rescue blanket Gold/Silver
- 2 Gauze 6cmx4m
- 3 Gauze 8cmx4m
- 1 Triangle bandage 96x96x136cm
- 1 Scissors
- 2 Wound wipes
- 4 Disposable Gloves vinyl
- 2 Disposable Medical Mask
- 1 First aid tutorial

Ref. No. 6500



Color: Green, Blue

EC - FIRST AID KITS

DIN Kit 13164 Medium



Table of contents

DIN 13164

- 1 piece of Tape 2,5cmx5m
- 14 Plaster assortment
- 1 First Aid Bandage 6x8cm
- 2 First Aid Bandages 8x10cm
- 1 First Aid Bandage 10x12cm
- 1 Burn compress 60x80cm
- 6 Compresses 10x10cm
- 1 Rescue blanket Gold/Silver
- 2 Gauze 6cmx4m
- 3 Gauze 8cmx4m
- 1 Triangle bandage 96x96x136cm
- 1 Scissors
- 2 Wound wipes
- 4 Disposable Gloves vinyl
- 2 Disposable Medical Mask
- 1 First aid tutorial

Other

- 1 Wound plastic spray 60 ml
- 1 Face shield with one way valve
- 2 Suture tape envelope 6,4x76mm
6 strips/envelope

Accessories: Wall bracket



Ref. No. 6120 + color

EC - FIRST AID KITS

Small



Ref. No. 6110 + color

Table of contents

- 1 Textile refill (36 plasters 25x72mm)
- 1 Textile Blue refill (36 plasters 25x72mm)
- 1 Textile refill (18 patches fingertip)
- 5 Wound wipes
- 1 First Aid Bandage 8x10cm

Accessories: Wall bracket

BURNCARE
RELIEVE BURN GEL

BURN PRODUCTS



Burncare 3,5 gr gel.

Tear open the bag and squeeze the contents over the burn. Let it be as thick a layer as possible. Don't rub it out. In case of fingertip injury: Put your finger in the pouch and leave it alone until the pain subsides. If necessary, seek medical advice.

Ref. No. BC5003

Burncare 50 ml and 120 ml gelspray.

Spray plenty of the contents over the burn. Let it be as thick a layer as possible. Do not rub out the gel. Cover the burn area with a wound compress. If necessary, seek medical advice.

Ref. No. BC5005 and BC50120



Burncare 5x5cm, 10x10cm, 20x20cm, 30x40cm (face mask) and 40x60cm compress

Tear open the bag and take out the compress and place over the burn. Squeeze out all the gel. Make sure that the compress covers the entire burn-area and 1-2 cm beyond this. Place the empty bag over the compress and fix with e.g. First aid bandage / alt. **SuperSoft**. Leave the compress on for 2 hours. **In case of injury to e.g. a hand: Put your hand in the bag and close the opening. In case of burns on the arm or leg: Open the package at both ends and thread through the arm or leg. Close the bag at both ends.** Larger/deeper burns, seek medical attention.

Ref. No. BC0505; BC5010; BC5020; BC5030 and BC5040

Burncare 25 gr gel tub

Squeeze the contents over the burn. Let it be as thick a layer as possible. Don't rub it out. Cover the burn area with a wound compress. If necessary, seek medical advice.

Ref. No. BC5025



BURNCARE BURN PRODUCTS

Burn Gel and Compresses
Can be customized as desired!

- Cools the burn
- Provides immediate pain relief
- Reduces the risk of infection



SMALL Ref. No. 5110

Content:
2 pcs 10x10cm compress
1pc Burncare spray 50ml
5 pcs 3.5gr gel sachets
1st First Aid Bandage

Accessories: wall bracket



MEDIUM Ref. No. 5120

Content:
1 20x20cm compress
2 pcs 10x10cm compress
1pc Burncare spray 50ml
5 pcs 3.5gr gel sachets
1st First Aid Bandage
2 shelves
1pcs wall bracket



10x10cm compress



BURNCARE BURN PRODUCTS

Burn Gel and Compresses
Can be customized as desired!

- Cools the burn
- Provides immediate pain relief
- Reduces the risk of infection



MEDIUM + FACEMASK
Ref. No. 5130

- Content:**
- 1 Facemask 40x30cm
 - 1 20x20cm compress
 - 2 10x10cm compress
 - 1 Burncare spray 50ml
 - 5 3.5gr gel sachets
 - 1 First Aid Bandage
 - 2 Shelves
 - 1 Wall bracket



10x10cm compress

MEDIUM - L Ref. No. 5070

- Content:**
- 1 Facemask 40x30cm
 - 2 20x20cm compress
 - 2 10x10cm compress
 - 1 Burncare spray 50ml
 - 5 3.5gr gel sachets
 - 1 First Aid Bandage
 - 2 Shelves
 - 1 Wall bracket



BURN - REFILL

Burn Gel- and compresses

- Cools the burn
- Provides immediate pain relief
- Reduces the risk of infection



Refill Small Ref. No. 5060

Content:

- 2 10x10cm compress
- 1 Burncare spray 50ml
- 5 3.5gr gel sachets
- 1 First Aid Bandage



Refill Medium Ref. No. 5061

Content:

- 1 20x20cm compress
- 2 10x10cm compress
- 1 Burncare spray 50ml
- 5 3.5gr gel sachets
- 1 First Aid Bandage



10x10cm compress



PHARMAPORE-PU

Waterproof plaster – post-operative plaster

POSTOPERATIVE WOUNDS - CUTS - LACERATIONS - ABRASIONS

Transparent protective plaster for mild to moderate secretion damage. Made from a polyurethane-, bacteria- and water-resistant film coated with adhesive sterilized with hypoallergenic Ethylene Oxide gas, supplied with non-adhesive centrally located wound pad.

- Large wound pad with non-adhesive coating
- Sterile
- Bacterial and water resistant
- Creates a humid atmosphere
- Provides very good wound healing
- Comfortable to wear / increases comfort
- Waterproof - allows for showering
- Allows free passage of gases and water vapour
- Self-adhesive polyurethane film around the wound pad

Clean the wound and remove the plaster from the packaging and Assemble according to the numbering on packages 1, 2 and 3.






*)

Sizes

Ref. No. 4267 8,25x6cm 200pcs/box

Ref. No. 4250 *) 7,2x5cm 200pcs/box **NEW**



Pharmapore – Sterile **NEW** Nonwoven Wound Dressings



Description

Pharmapore sterile combines a low-adhesion, highly absorbent wound pad with a soft malleable material for an easy and effective treatment of wounds.

Function

- Individually packed in a peel-off bag and sterilized with EO gas.
- Soft, breathable dressing.
- Highly absorbent and malleable wound pad.
- Hypoallergenic glue.

Benefits

- Highly absorbent, low adhesive pad.
- The polyethylene wound contact layer allows the passage of blood or secretions to a highly absorbent wound pad.
- While maintaining a good healing environment and minimal pain during removal.

Safe and secure fixation

- Hypoallergenic adhesive that allows the skin to breathe to protect the skin from maceration.
- This ensures a safe and secure fixation of the dressing.

Soft and malleable

- Soft and conformable nonwoven material ensures patient comfort and easy application to hard-to-reach and non-smooth wound sites, such as forearm and shoulders, hands and fingers

Patient Comfort

- Thick wound pad prevents fluid from penetrating and protects the wound from trauma.

Use

- Pharmapore sterile is ideal for post-operative use, cuts, lacerations and sutured wounds.



7.2 x 5cm 200pcs / box
Ref. No. nr: 4350

8.3 x 6cm 200pcs / box
Ref. No: 4360



STERILE EO



PHARMASTRIP


Suture Tape – Cuts

Att! 6pcs strips/envelope

Size: 6,4x76mm

Flexible skin sealing strips made of polyamide, non-elastic backing coated with a pressure-sensitive, hypoallergenic adhesive.

- Breathable, porous, non-elastic carrier can be used on contoured areas on the body or where bending and movement occurs.
- White color for a cosmetically pleasing look.
- Eliminates marks such as stitches and staples,
- A better cosmetic result than sutures or staples.
- Quick and easy to apply to save time.
- Hypoallergenic adhesive.
- Packed in sealed bags E.O. sterilized.
- For use in highly profiled areas or areas of the body with movements of muscles and joints.
- Provides support to wounds after early suture or staple removal.
- Useful for supporting wounds postoperatively and after suture removal.



Sizes

Ref. No. 6476-25 25pcs envelope / bag
Ref. No. 6476 100pcs envelope / box

STERILE EO


CE
0459



ELASTIC ADHESIVE BANDAGE

7,5cm x 5m

- Longitudinal Elastic Adhesive Bandage
- High quality 100% stretchable cotton textile
- Porous bandages
- Provides a stable but controlled support
- Porous bandage allows sweat to evaporate to avoid skin maceration
- Soft edges prevent the bandage from creasing or leaving marks on the skin
- Retains elasticity
- Varicose veins inflammation and ulcers (phlebitic and thrombotic)
- Sprains & strains
- Sports Injuries
- Prevention, treatment and rehabilitation of injuries



Ref. No. 4750 7,5cm x 4,5m

CE



RESCUE BLANKET

Gold / Silver

160x210cm



The life-saving blanket can be used, for example, in industry, in traffic, workplaces and leisure. As well as very suitable even in case of accidents, burns and bad weather. However, the face should not be covered.



RÄDDNINGSFILT GOLD / SILVER

Den livräddande filten kan t.ex. användas inom industri, trafiken, arbetsplatser och fritid. Vid olyckor, brännskador och vid oväder. Ansiktet skall inte täckas över.



160 x 210 cm

Skydd mot kyla: Silverfärgade sidan måste vändas inåt. Den Aluminiserade ytan reflekterar uppåt 80% av kroppstemperaturen. Hypotermi och frystning förebyggs på detta sätt.

Skydd mot värme: Silverfärgade ytan måste vändas utåt. Solstrålar reflekteras och håller kroppen sval (eller till bilen, tält, etc.)



Use

Protection from the cold: The silver-colored side should be turned inwards. The silver-colored surface reflects up to 80% of the body temperature. Hypothermia and freezing are prevented in this way.

Protection against heat: The silvery surface should be turned outwards. The sun's rays are reflected and keep the body cool. Can also be used for the car, tent, etc.

Ref. No. 3161

160x210cm



BLOODSTOPPING - BANDAGE

Gauze / Wound compress

Blood stopping bandage is a first aid kit. It is a combination of a silicone-treated wound pad to avoid getting stuck in wounds, and a porous gauze bandage of 400cm. A practical dressing for the treatment of bleeding wounds.

Blood stopping bandage can also be used as a pressure bandage for major bleeding by combining with an extra gauze bandage over the bleeding.

Sizes

Ref. No. 3158 10x12 cm

Ref. No. 3159 8x10 cm



STERILE | EO



VENDAPRESS Self-adhesive flexible bandage

Elastic bandage to support sprains

Product

This product is constructed from a combination of the elastic and high strength of the polyester fiber and a latex coating that allows for surfaces that are self-adhesive. You can't feel any stickiness on either side of the bandage, but it sticks completely to itself. Does not stick to skin, fur or wounds.

Benefits

Tear off, it can be torn off very easily and smoothly

Can be used in different thicknesses and for all body parts if needed

Self-adhesive, no need for staples or tape

Retains its original size – does not contract

Elasticity, very good elasticity, and good elasticity

The bandage is breathable

Quick and easy to remove, leaves no adhesive residue

Use

Fixation of e.g. burn compresses

Elastic dressing, plaster dressing and protective dressing during assembly work.

First aid product for sprains and sprains



Ref. No.	Width	Length
8133-blue	7,5 cm	4,5 m

LONG PLASTERS

Long plasters– 2x18cm x 100pcs


Long Plaster is a wound dressing in elastic textile with absorbent wound pad. The long plaster means that, in addition to covering the wound with the wound pad, you can also cover scratches next to the wound. The patch is available as Blue Detectable Textile and as Textile.

- Highly absorbent but with low adhesion wound pad
- Highly hypoallergenic acrylic adhesive for skin sensitivity
- Elastic and comfortable to use
- Painless removal of the patch

Sizes 100pcs/box

Ref. No. 3361-18/2 Textile - elastic

Ref. No. 3360-18/2 Blue Detectable Textile - elastic



Sensible Tape

with dispenser - 2,5cm x 4,55m



Description

Tear-off paper tape in a material with low allergy risk. Self-adhesive tape made of micro-porous non-woven material coated on one side with a hypoallergenic, skin-friendly adhesive.

Properties

- Non-woven material.
- Skin-contact adhesive
- Workable
- Latex-free
- Less pain during removal
- Permeable to air and water vapour

Use

- Fixation of bandages and compresses
- Use for sensitive and dry skin where ordinary tape cannot be used
- Apply to clean and dry skin

Storage

- Keep away from light and heat
- Do not store above 25°C

Packaging	6pcs / box
Ref. No.	32545

DISPENSER - WOUND WIPE

Contains 45 wound wipes



- Transparent Plastic Dispenser
- Contains 45 wipes for wound washing
- Wall mounting for easy use
- Easy to see how many napkins are remains in the dispenser
- Easy refilling of napkins

Ref. No. 3201-D

STERILE EO

CE
0459



Wound wipe

7,5x14cm 45pcs



- Alcohol free Wound Wipes - Refill
- Size 7,5x14cm
- Contains 1% Cetrimid which is an effective antiseptic against bacteria and fungi
- Tear open the edge of the package and remove the Napkin

Ref. No. 3200

STERILE EO

CE
0459



WALL BRACKET - SPRAY

Wall bracket for sprays



Wall bracket – small 35mm

Suitable for spray bottles with a diameter of 35mm for example.
Material – PET plastic


- Blood Stop Spray, 60ml
- Wound Wrap Spray, 60ml
- Burncare spray, 50ml
- Eye Wash Spray 50ml

Ref. No. 2040

Wall bracket – large 52mm

Suitable for spray bottles with a diameter of 52mm e.g. Eyewash spray 250ml
Material – PET plastic

Ref. No. 2050



Triple Holder – 52mm

Suitable for spray bottles with a diameter of 52mm e.g. Eye wash spray

Material – PET plastic

Ref. No. 2060

FIRST AID KITS

With wall mount

Size: 31x21x12cm

A practical first aid bag in impact-resistant ABS hard plastic with wall mount. This allows one to quickly take the bag to an injured person for a quick treatment.

The case has double inner lids to keep the products in place. Thanks to its clarity, it is easy to see what is available and when it needs to be replenished.

Design the case yourself for your own workspace by choosing products outside the standard range or only with products outside Din 13164.

Suitable for both large and small workplaces, sports arenas, boats and ships, the transport industry, rescue vehicles, etc.

Table of contents

DIN 13164

- 1 Piece of Tape 2,5cmx5m
- 14 Plaster assortment
- 1 First Aid Bandage 6x8cm
- 2 First Aid Bandages 8x10cm
- 1 First Aid Bandage 10x12cm
- 1 Burn compress 60x80cm
- 6 Compresses 10x10cm
- 1 Rescue blanket Gold/Silver
- 2 Gauze 6cmx4m
- 3 Gauze 8cmx4m
- 1 Triangle bandage 96x96x136cm
- 1 Scissors
- 2 Wound wipes
- 4 Disposable Gloves vinyl
- 2 Disposable Medical Mask
- 1 First aid tutorial



Ref. No. 3165



Ref. No. 3160 Complete
3170-G Empty
Color: Green

FIRST AID KIT

DIN – KIT 13164

43 delar

A first aid kit with contents for the treatment of wounds. Suitable for most workplaces.

DIN Kit 13164 is a very suitable standard range that can be placed in AlphaMedic bags or in existing First Aid kits, boxes and pillows.

Suitable for both large and small workplaces, sports arenas, boats and ships, the transport industry, rescue vehicles, etc.

Ref. No. 3165

Table of contents

DIN 13164

- 1 Piece of Tape 2,5cmx5m
- 14 Plaster assortment
- 1 First Aid Bandage 6x8cm
- 2 First Aid Bandages 8x10cm
- 1 First Aid Bandage 10x12cm
- 1 Burn compress 60x80cm
- 6 Compresses 10x10cm
- 1 Rescue blanket Gold/Silver
- 2 Gauze 6cmx4m
- 3 Gauze 8cmx4m
- 1 Triangle bandage 96x96x136cm
- 1 Scissors
- 2 Wound wipes
- 4 Disposable Gloves vinyl
- 2 Disposable Medical Mask
- 1 First aid tutorial



Feel free to contact us:

- ❖ If you want to place an order. See addresses below!
- ❖ If you have products that need to be customized to best suit the market and your customers!
- ❖ If you want more info about products and retailer prices!
- ❖ If you have ideas that you want to bounce ideas off us!

Yours sincerely,
AlphraMedic AB



Addr: Kämpevägen 17 - 553 02 Jönköping

Phone: + 46 36 340485

Mob: +46 70 273 77 40

Jan

Mob: +46 70 993 99 10

Patrik

E-mail: info@alphramedic.se

Web: www.alphramedic.se