


Pharmagel Comfort

**CHRONIC WOUNDS - SMALL BURN -
SORE - PRESSURE ULCERS –
TREATMENT AFTER THE BURN**



Description:

Pharmagel Comfort is a sterile transparent polyurethane gel dressing containing 60% water, with hypoallergenic adhesive polyurethane Film border to prevent bacterial penetration.

Benefits:




- Creates a moist environment, suitable for treating chronic wounds, promoting wound healing, while slight cooling effect aids patients comfort..
- Allows free passages of gases and moisture vapour, protecting skin against maceration.
- Leaves no residue on wound upon removal of the dressing.
- Pharmagel is transparent, allows healing process to be observed without removing the dressing, promoting formation of new tissue during granulation phase.
- Pharmagel doesn't dissolve by absorbed wound exudates.
- Pharmagel Comfort are individually packed in moisture proof pouches sterilized by radiation.

Indications:

- Suitable for treatment of chronic and slow healing wounds in granulation epithelialisation phase.
- Due the three dimensional structure of Pharmagel and high water content it is suitable for use in debriding and for continuous treatment on second degree burns after using Burnaid gel.

Contraindications:

Pharmagel should not be used on infected wounds nor on third degree burns, it should not be used also on highly exudation wounds nor on exposed muscle, bones nor tendons.



<u>Ref. no:</u>	<u>Size:</u>
4565-C	45x65mm
4575-C	75x100mm

