

# PLASTIC DRESSING -spray



CE  
2265

60ml 80



## Liquid plastic


- ◆ a plastic wound dressing based on acryl copolymer for an easy and quick treatment of lacerations, scratches and minor injuries
- ◆ after application, it dries within 3 minutes
- ◆ it creates a waterproof protecting film
- ◆ the dressing is not water-soluble
- ◆ it stops micro- organism from surrounding environment
- ◆ *it does not stop natural exudates*
- ◆ the wound remain closed against disadvantaged effect, still the healing process continues under the film dressing
- ◆ the dressing disappear by it self after 3-4 days by gradually abrasion of the plastic film

## Protective film for the treatment of surface wounds

- providing protection for minor wounds such as lacerations, abrasions and cuts
- reduces the risk of infection, protects against water and dust
- good in areas where it is difficult to use ordinary plaster
- spray a thin layer, about 15-20cm from the wound, repeat if necessary

Contains: Ethanolum denaturatum, Dimethyletherum, Copolymerum acrylatum, Dioxacyclopentanum.

Art. Nr: 3140



AlphaMedic