





PRODUCTS



DISPENSER

to 5 meters rolls

Suitable for all types of 5m patch rolls

D-56 is a popular dispenser for a simple way to cut the size of the patch one needs, wide or narrow, your choice!

A Practical, Easy, Economical and flexible system! *Mounted on the wall!*

Made of clear acrylic plastic. The dispenser are available in three different widths: 3cm, 6cm and 8cm.

Ref. no. 3000-3 D-53 to 3 cm x 5m

Ref. no. 3000-1 D-56 to 6 cm x 5m

Ref. no. 3000-8 D-58 to 8 cm x 5m



Suitable for all types of plaster rolls ex. Textile, Non Woven, Polyurethane Foam, etc.

ADHESIVE PLASTER ON ROLL

Wide: 6cm and 8cm

Textile Plasters

- Use on clean and dry skin
- Elastic and comfortable to use
- Painless removal of the patch
- Fits D-56 patch dispenser
- Cut to desired size

Ref. no. 3361 6cm



A powerful Textile Plasters

- Use on clean and dry skin
- Elastic and comfortable to use
- Painless removal of the patch
- Fits D-58 patch dispenser
- Cut to desired size

Ref. no. 3361-8 8cm

Non Woven patches, Hypoallergenic

- Suitable for sensitive skin
- Use on clean and dry skin
- Elastic and comfortable to use
- Painless removal of the patch
- Fits D-56 or D-58 patch dispenser
- Cut to desired size

Ref. no 3362 6cm; 3362-8 8cm



ADHESIVE PLASTER ON ROLL

Wide: 6cm and 8cm

Blue Detectable Textiles with metal foil

- Use on clean and dry skin
- The blue color provides visual identification
- wound pad with metal foil can be identified by metal detector
- Elastic and comfortable to use
- Painless removal of the patch
- Fits D-56 or D-58 patch dispenser
- Cut to desired size

Ref. no 3360 6cm; 3360-8 8cm



Foam patch

- Micro perforated
- Use on clean and dry skin
- Elastic and comfortable to use
- Painless removal of the patch
- Fits D-56 or D-58 patch dispenser
- Cut to desired size

Ref. No 3364 6cm

SUPERSOFT



6cm x 5m

SuperSoft

Blue

- Self-adhesive foam dressing
- The blue color provides visual identification
- Elastic and comfortable to use
- Does not stick to the wound or in the hair
- The dressing is breathable and provides great healing
- Compatible with D-56 patch dispenser
- Cut to desired size
- For fixation of burn dressings

Art. nr 3385-B (blue)



SuperSoft

Neutral

- Self-adhesive foam dressing
- The blue color provides visual identification
- Elastic and comfortable to use
- Does not stick to the wound or in the hair
- Compatible with D-56 patch dispenser
- Cut to desired size
- For fixation of burn dressings
- The dressing is breathable and provides great healing

Art. nr 3385 (Neutral)



EASYON DISPENSER

Plaster dispenser with plasters refill



Ref no.2505 - complete

EasyOn plaster are sterile, easy and comfortable to use!!



Ref no.2510- empty

EasyOn dispenser is made of transparent and durable plastic. To be used with EasyOn plasters refill. See below. With a simple operation, replace the refill. Mounted on the wall.

EasyOn refill fits EasyOn dispenser and a single action, replace the refill. Gently pull down the lower edge of the patch to loosen the top of the paper on the patch back - ready to assemble!

EasyOn is supplied in three different material and three different sizes.

PU plastic – elastic Blue Detectable Fabric Elastic fabric



Ref no.	Product	Size	Plaster/refill
3358-25-1R	PU-plast Blue Detec	25x72mm	40 New
3359-25-1R	PU-plast	25x72mm	40
3359-52-1R	"	52x72mm	20
3360-25-1R	Blue Detec Textil	25x72mm	36
3360-52-1R	"	52x72mm	18
3360-FT-1R	" "	Fingertopp	18
3361-25-1R	Textilplåster	25x72mm	36
3361-52-1R	"	52x72mm	18
3361-FT-1R	"	Fingertopp	18



EasyOn patches are sterile, simple and convenient to use. Opportunity to own refill are at larger quantities!



PLASTIC DRESSING

Plastic dressing spray is an aerosol with quick drying wounds plastic that provides a fast and reliable protection. It is designed to protect minor wounds and injuries such as lacerations, abrasions and lacerations on sports, work-places and homes.


Plastic dressing spray reduces the risk of infection, protect against water and dirt which allows safe showering and stay in wet environments or dusty and dirty work, easy, convenient and safe to use, great in areas where it is difficult to use traditional plasters and dressings.

Ref no. 3140

80ML



CE
0120



NEW

BLOODSTOP SPRAY alginate

BLOODSTOP spray is an aerosol with blood-stopping powder. It is designed to stop minor bleeding from superficial injuries such as lacerations, abrasions and lacerations on sports, workplaces and homes.

Press the spray button to apply the haemostatic powder to the wound tissue that mainly gets disinfected. Spray from a distance of about 15-20 cm. Press the spray button for about 2 seconds, by repetitive treatment shake the bottle thoroughly each time.

CE
1023

Ref no. 3130

80ML




PHARMAGEL

Moist wound healing



CHRONIC WOUNDS - SMALL BURN - SORE - PRESSURE ULCERS – TREATMENT AFTER THE BURN

- The gel pad contains 60% water
- Cools and soothes
- Creates a moist atmosphere
- Provides very good wound healing
- Minimizes scarring and sloughing
- Accelerates the formation of new tissue
- Comfortable cushioning against blisters / sores
- Comfortable to wear / increase comfort
- Waterproof - allows showering
- Allows free passage of gases and water vapor
- Self-adhesive polyurethane film around the gel pad



treatment of pressure ulcers, not on exuding wounds!

Sizes

Ref. no 4565-C 45x65mm

Ref. no 4575-C 75x100mm

BAGG FIRST AID BAGS



- Initial design / Registered Trademark
- Durable, lightweight 420 denier polyester.
- Weather and water resistant material.
- Optional mounting points (daisy chain).
- Double Sided "Heavy Duty" zipper for easy opening and closing.
- Available in 4 different eye-catching colors: red, yellow, orange and green.
- Robust plastic handle with extra grip.
- Reinforced belt anchor.
- Extra "Molle" attachment point.
- Interior divided into several sections to keep the contents in place.
2 large sealable sections with zipper and two elastic sections.
- Dimensions: 21 x 15 x 11cm.



Bagg first aid bags are supplied filled or empty. Content can easily be adapted to suit your customers' needs.

Art no. 5510 - color (empty)

No. 5515 - color (complete with DIN kit 13164)

EC – FIRST AID CASES

Wounds - Burns - Eye Injury

Stylish and practical quality bags that are supplied complete, customized or empty. Perfect in the workplace where you can build up a complete working system for first aid treatment. The bags can be mounted beside or under each other.

When wall mounted it takes a single action to get the whole case to the injured. Choose shelves for clarity. Can also be used without wall bracket as necessary.

The bags can also be used for other products such as gloves, protective goggles, dust mask hearing protection, etc..

SIZES

Large, Medium; Small

DIMENSIONS:

Large: 30x33x12

Medium: 30x23x12

Small: 15x23x8

ACCESSORIES:

Bracket L/M, S

Shelves L, M

COLOUR:

Green, Blue, Red (Burn Kit)

Transparent (Eyewash Kit)



EC – FIRST AID CASES

DIN+EASYON

Large and Medium



CONTENTS:

DIN KIT 13164 + EasyOn dispenser with 2 plasters refill incl. 2 shelves

Accessories: Wall Bracket

DIN KIT

EasyOn dispenser

1 Textile-refill + 1 PU-plastic refill
25x72mm 25x72mm

Table of Contents

DIN 13164

- 1 Tape 2.5 cmx5m
- 14 Plaster assortment
- 8 Non-woven patch 10p 6x50cm
- 3 Ambulance Dressing 8x10cm
- 1 Ambulance Dressing 10x12cm
- 2 Burn dressing 40x60cm
- 1 Burn dressing 60x80cm
- 6 Swabs 10x10cm
- 1 Rescue Blanket Gold / Silver
- 2 Gauze 6cmx4m
- 3 Gauze 8cmx4m
- 2 Triangular bandages 96x96x136cm
- 1 Scissors
- 2 Alcohol free wipe
- 4 Disposable gloves vinyl



Colour: Green, Blue

Burnaid®

BURN GEL PRODUCTS

Burn Gel- and Dressing
Can be adapted to your wishes!

- Cools the burn
- Provides immediate pain relief
- Reduces the risk of infection



Medium



Small



10x10cm



The nylon bag for larger dressings ex. (facemask) 40x30 cm and 55x40cm or KIT - L and KIT - XL

Burnaid®

BURNAID PRODUCTS

- Cools the burn
- Provides immediate relief
- Reduces the risk of infection



3.5g gel



10x10cm



20x20cm



120x200cm body



55x40cm



40x30cm facemask



50 ml spray

Dressings:

10x10, 20x20, 30x40 (face mask), 55x40 and 120x200cm (body)

Gel:

3.5g sachets, 50ml pump spray

PHARMAPORE-PU

Waterproof plasters - postoperative plasters

Postoperative wounds - Cut - Laceration - Graze




Clear protective plasters for light-to-moderate secretions damage. Made of polyurethane, bacterial and waterproof film coated with hypoallergenic adhesive sterilized with ethylene oxide gas, is supplied with a non-adherent wound pad located centrally.

- Large wound pad with non stick coating
- Sterile
- Bacterial and water safe
- Creates a moist atmosphere
- Gives a very good wound healing
- Comfortable to wear / increase comfort
- Waterproof - allows showering
- Allows free passage of gases and water vapor
- Self-adhesive polyurethane film around the wound pad

Sizes:

Ref. no 4767 8,25x6cm 200pcs/box

Ref. no 4712 8,25x12cm 60pcs/box



Clean the wound and remove the patch from its package and mount according to the numbering of the pack 1, 2 and 3

PHARMACOLL

Hydrocolloid plasters - blister plasters

CHAFING- PRESSURE ULCERS - BLISTERS - EXUDING WOUNDS

Pharmacoll is a blister plasters, designed according to modern and effective wound healing. Self stitched and elastic. A hydrocolloid dressing for blisters and chafing. They are individually wrapped and sterile.

Pharmacoll seems analgesic when hydrocolloid bearing particles form a gel cushion in absorbing moisture and wound secretion which provides a moist wound environment, which promotes the formation of new skin cells. Protect nerve ends to pressure and shocks.

Pharmacoll does not stick to the wound and the dressing can be easily and painlessly removed and without affecting healing and tissue.

- Sterile
- High absorbency
- Bacterial and water safe
- Creates a moist atmosphere
- Gives a very good wound healing
- Comfortable to wear / increase comfort
- Waterproof - allows showering
- Minimizes crust and scarring
- Allows free passage of gases and water vapor



*Clean the wound thoroughly. **Pharmacoll** can adhere properly only if the affected area is dry and free of lotions. Remove the backing and press firmly, put your hand over and heat for about 1 minute for the edges of the patch to adhere properly.*

***Pharmacoll** should sit on until it comes off by itself after a few days*

Sizes:

Ref. no 3382

68x43mm

50pcs/box

PHARMASTRIP


Suture tape - cuts

6x strips / envelope

6,4x76mm

Flexible skin closure strips made of polyamide fabric, non elastic backing coated with a pressure-sensitive, hypoallergenic adhesive.

- Breathable, porous, non-elastic backing, can be used on contoured body areas or where bending and movement occur.
- White or tan color for a cosmetically pleasing appearance.
- Eliminates crosshatching and puncture marks, gives a better cosmetic result than sutures or staples.
- Fast and easy to apply for to save time.
- Hypoallergenic adhesive.
- Packaged in sealed bags
- Available in sealed pouches sterilized.
- For use in highly contoured areas or body areas of musculoskeletal movements like joints.
- Providing wound support following early suture or staple removal.
- Useful for wound support post operatively and following suture removal.



Sizes

Ref. no 6476-25 6,4x76mm 25pcs envelope / zipbag

Ref. no 6476 6,4x76mm 100pcs envelope / box



CE
0120

STERILE EO



RESCUE BLANKET

Gold / Silver

160x210cm



The life-saving blanket can example. be used in industry, in transport, workplaces and leisure. And highly appropriate also to accidents, burns, and storms. The face should not be covered.



Use

Cold protection: The silver surface should face inwards. The silver surface reflects upwards of 80% of body temperature. Hypothermia and freezing is prevented in this way.

Protection against heat: The silver surface should face outward. The sun's rays are reflected and keeps the body cool. Can conveniently be used to the car, tent, etc.

Ref. no 3161 160x210cm

AMBULANCE DRESSING

Gauze / Wound dressing

Ambulance dressing is a first aid dressing. It is a combination of a silicone-treated wound pad to avoid getting stuck in the wound, and a porous gauze 400cm. A practical dressings for treatment of bleeding wounds.

Ambulance dressing can also be used as a tourniquet for major bleeding by combining with an extra gauze over the bleeding.

Sizes

Ref. no 3158 10x12 cm

Ref. no 3159 8x10 cm



PHARMAFIX - PU


Polyurethane film - protective film on roll

Pharmafix - PU is a transparent adhesive film on roll, made of polyurethane with an adhesive of acrylic. It is specially designed for dressing retention and protection of wounds, such as strengthening the joint.

Water and bacteria proof

The film is impervious to water and acts as a bacterial barrier which means that the patient can shower and swim with the dressing and to be staying in humid or wet environments. The permeability to water vapour causes the excess moisture to evaporate out. The skin-friendly adhesive surfaces of acrylic reduces the risk of irritation when removing the film.

Ref. no 4100 5cmx10m; Ref. no 4125 5cmx2,5m



2,5m and 10m roll



LONG ADHESIVE PLASTER

Long patch elastic – 2x18cm x 100pcs


Long plasters is a wound dressing in elastic fabric with absorbent wound pad. The long patch allows one to in addition to covering the wound with wound pad also cover scratches next to the wound. The patch is available as Blue Detectable Textiles and as Textile.

- Superabsorbent but with low adhesive wound pad
- Very hypoallergenic acrylic adhesive skin sensitivity
- Flexible and convenient to use
- Painless removal of the patch

Sizes 100pcs/box

Ref. no 3360-18/2 Blue Detectable Textile

Ref. no 3361-18/2 Textile



SURGICAL TAPE – TEXTILE

Fixation of all types of dressings



- Strong Surgical Tape for normal skin.
- Coated with zinc oxide rubber adhesive with high initial and permanent adhesion at where it attaches strong even under heavy stress.
- High tensile strength.
- Zigzag edges making it easy to tear by hand.
- The product has been made permeable to air and moisture through a special process which avoids the possibility of maceration due to the exclusion of air.

Ref. no 3361-KT25

2,5 cm x 4,55 m

Vendapress

Elastic and self-adhesive
7,5cm x 4,5m

Advantages

- Tearing off, it can be torn off very easily and smoothly
- Self-adhesive, no need for staples or tape
- Retains its original size - not deducted together
- Elasticity, very good elasticity and good ductility
- The bandage breathe
- Quick and easy to remove, leaving no sticky residue

Ref. no 8133-blue

Product

This product is constructed from a combination of polyester fibers stretch and high strength and a latex coating that allows surfaces that are self-adhesive. One can not experience any stickiness on any side of the bandage, but it adheres completely to it self.



DISPENSER – WIPES

Contains 35pcs - alcohol free wipes



- Transparent plastic dispenser
- Includes 35 wipes for cleaning wounds
- Wall mounting for ease of use
- Easy to see the number of napkins remaining in the dispenser
- Easy refill of wipes

Ref. no 3201-D

ALCOHOL FREE WIPE

10x14cm 35pcs



- Alcohol free wipe- Refill
- Size 10x14cm
- Contains 1% Cetrimide which is an effective antiseptic to bacteria and fungi
- Tear up the edge of the box and take out the napkin

Ref. no 3200



WALL BRACKET – SPRAY

Wall bracket for sprays



Wall bracket - Small 35mm

Suitable for spray bottles with 35mm diameter ex. Blood Stop spray, Wound plastic spray and Burnaid spray

Material - PET plastic

Ref. no 2040

Wall bracket – large 52mm

Suitable for spray bottles with 52mm diameter ex. Eyewash spray

Material – PET plastic

Ref. no 2050



Wall bracket – mushroom cap

Fits eyewash spray bottles with mushroom cap - triple. Acrylic

In use - pull the bottle straight down so the cap comes off and remains in the wall bracket. Space for 3x bottles.

Ref. no 3220

